Verso

User Manual

CONTENTS

ENGLISH

USER MANUAL 3

ESPAÑOL

MANUAL DE INSTRUCCIONES 21

PORTUGUÊS

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FRANÇAIS

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DEUTSCH

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DUTCH

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ITALIANO

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MODEL NO.: 9P-046540

PLEASE READ THE MANUAL BEFORE USE

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USING

DANGER – To reduce the risk of electrocution:

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store this product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately at the mains power point.

WARNING – To reduce the risk of burns, electrocution, fire, or injury to persons:

1. Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
2. This system is not for use with patients who have a spinal cord injury.
3. Close supervision is necessary when this product is used for or near children. Electrical burns or choking incident may result from a child swallowing a small part detached from the device.
4. Use this product only for its intended use as described in this manual. Do not use other mattresses not recommended by the manufacturer.
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your supplier or Apex Medical Corp. for examination and repair.
6. Keep the cord away from heated surfaces.
7. Never block any air openings of this product or place it on soft surfaces, such as a bed or couch, where openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
8. Never drop or insert any object into any opening or hose.
9. Do not modify this equipment without authorization of the manufacturer.
10. Mattress covers have passed skin sensitization and skin irritation test. However, If you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
11. Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.
12. Keep the pump away from flammable liquids or gases.
13. Do not use the extension cord for a long time
14. Thrid conductor in the POWER SUPPLY CORD is only a functional earth.
15. The Operator cannot operate the instrument and touch the patient at the same time.

CAUTION –

If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

NOTE, CAUTION AND WARNING STATEMENTS:

NOTE – Indicates useful tips.

CAUTION – Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING –Raises awareness to a potential danger that requires correct procedures or practices in order to prevent personal injury.

SYMBOLS

	Authorized representative in the European community.
	Catalog, reorder or reference number
	Manufacturer
	Complies with standards protecting against electric shock for type BF equipment.
	Attention, you should read the accompanying information carefully!
	Consult operating instructions for use
IP21 	Protected against solid foreign objects of 12.5 mm and greater; Protection against vertically falling water drops
	Class II
	Temperature limitation/temperature range
	Dry clean, Any Solvent Except Trichloroethylene
	Do Not Iron
	Tumble Dry, Normal, Low Heat
	Do Not Tumble Dry
	Do Not Bleach
	Do Not Dry Clean
	Machine wash, regular / normal, 95 degrees C (203 degrees F)
	Machine wash, regular / normal, 60 degrees C (140 degrees F)
	Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.
	CE Marketing
	"ON" (power)
	"OFF"(power)

1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

1.1 General Information

The system is a high quality and affordable mattress system suitable for treatment and prevention of pressure injuries.

The system has been tested and successfully approved to the following standards:

IEC/EN 60601-1
IEC/EN 60601-1-2
IEC/EN 61000-3-2 Class A
IEC/EN 61000-3-3
CISPR 11 Group 1, Class B

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

1.2 Intended Use

This product is intended to help and reduce the incidence of pressure injuries while optimizing patient's comfort. It can support patient up to 160 kg (5" overlay) and 180 kg (8" replacement). And It can also provide the following purposes:

- For long term care patients who suffering from pressure injuries.
- For pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in the knowledge of prevention and treatment of pressure injuries.


NOTE: Equipment is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

2. Product Description

2.1 Pump and Mattress System


1. CPR
2. Pump Unit
3. Mattress System
4. Quick Connector
5. Transport Cap

2.2 Pump Unit


Front

1. Quick Connector Slot
2. Display Panel
3. Transport Cap


Rear

4. Power Cord
5. Hangers
6. Fuses
7. Filter

2.3 Display panel

1. Pressure Adjustment Knob

Pressure adjustment knob controls the air pressure output. When turning clockwise, the output pressure will be increased. On the other hand, when turning counter-clockwise, the output air pressure will be decreased. The pressure adjustment indicator would represent different weight support; however, it mainly depends on the type of mattress.

2. Low Pressure Indicator

When low pressure LED lights up, the pressure inside of air mattress is below normal. Please refer to troubleshooting.

3. Normal Pressure Indicator

When green LED lights up, the pressure inside of air mattress green reaches the preset pressure setting.


4. Main Power Switch


To turn on/off the pump unit with illuminated green light.

5. Therapy Mode Button


Switch between Static and Alternating mode of mattress.

2.4 Weight Indicator Of Verso 5" And 8" Systems:

2.5 Lengthened hangers (Optional)

1. Lengthen Pump hanger (Optional)

We, APEX MEDICAL, also offer a lengthened hangers to fit different types of nursing bed.


3. Installation

Unpack the box to check if there is any damage which may have occurred during shipment. If there is anything damaged or missing, please contact your local service center immediately.

Package Content List

- Mattress unit x 1 (may not be included if only pump unit is purchased)
- Pump unit x 1
- User manual x 1

3.1 Pump & Mattress Installation

1. Place the mattress or pad on top of a bed frame. Orientate the mattress so that the foot-markers on the mattress are situated towards the end of the bed (where pump will be located).
2. Place the pump hangers over the bed railing at the foot of the bed. Support the pump from the bottom and turn the hangers outward (or inward) to secure the pump against the railing.
3. Connect quick connector from air mattress to the pump unit. When a "click" is felt or heard, the connection is completed and secured.

⚠ NOTE: Check and ensure the air hoses are not kinked or tucked under mattress.

4. Plug the power cord into electrical outlet.

⚠ NOTE: Make sure the pump unit is suitable for the local power voltage.

5. Turn the main power switch on the display panel to ON  position.

⚠ NOTE: You can unplug the unit to power-off the device.

⚠ CAUTION: The pump supplied with the mattress can only be used for mattresses recommended by the manufacturer. Do not use it for any other purpose. (applied part: air mattress)


⚠ NOTE: During power outages, you can cover the Quick Connector with the transport cap to maintain air pressure inside the air cells.

⚠ NOTE: After installation, make sure any excess cord or tubing is removed from possible foot traffic to avoid accidental tripping. All EQUIPMENT should be positioned to always allow unhindered patient access by physicians and caregivers.

The EQUIPMENT should be firmly placed at position where users/ caregivers can access easily.

4. Operation

4.1 General operation


1. Turn the main power switch to "⏻".
2. Once the pump begins delivering air into the mattress it will take approximately 30 minutes to fully inflate the mattress.
3. The pump unit is beginning to deliver air into the mattress. Turn the pressure adjustment knob clock-wise to increase firmness (inner pressure). When the initial inflation is completed, the system will enter the therapy mode that was chosen by users.
4. The low pressure indicator (yellow LED) will light up when the mattress is not yet fully inflated.
5. When the appropriate pressure is reached, the low pressure indicator (yellow LED) will go off and normal pressure indicator (green LED) will light up automatically.

NOTE: Every time when mattress is first setup for use, it is recommended that the pressure adjustment knob be set to "Firm" for the quickest inflation. User can then adjust air mattress to the desired softness afterward.

4.2 Emergency CPR Operations

CPR must be performed on a firm surface. Therefore, if an emergency CPR situation occurs with the patient on the mattress, the mattress must be quickly deflated. This is done by quickly pulling on the CPR tag located at the head of the mattress on the patient's right-side. The quick connector found from the pump unit can also be disconnected to hasten deflation.

4.3 Pressure Set Up

Users can adjust the pressure of air mattress to a desired support by adjusting the pressure adjustment knob. Please consult your physician for a suitable setting. The pressure indicator would represent different weight support; however, it mainly depends on the type of mattress.

A simple guideline of mattress type

Mattress Type	5" mattress	8" mattress
Max. Support	160 Kg	180 Kg

When a Normal Pressure indicator lights up, pressure is reached to the desired firmness that patient can lie on the mattress.


⚠ NOTE: Check to see if the suitable pressure is selected by sliding one hand between the deflated air cells and the patient to ensure patient is not bottoming out. Users should be able to feel the minimum contact.

4.4 Low Pressure Indicator

When an abnormal low pressure occurs, the Low Pressure indicator (yellow LED) will light up. Check that connections are correctly made and that they are correctly installed as per installation instructions.


⚠ NOTE: If the pressure level is consistently low, check for any leakage (tubes or connecting hoses). If necessary, replace any damaged tubes or hoses. Or contact local qualified service center (dealer) for repair.

4.5 Therapy Mode Selection

Users can press the Therapy Mode button to select between Alternating mode and Static mode. In the Alternating mode, the air cells are partially deflated and inflated periodically to avoid prolonged pressure on any single point of the patient; this is to prevent pressure injuries from happening. In the Static mode, all of the air cells are equally inflated to provide immersion and envelopment to the patient.

4.6 Transport Mode

When there is a power outage or a need to transport the patient, the quick connector can be capped as shown to slow down the air leakage. Depending on the setting and mattress height, the mattress can remain inflated for a period of time. Verso 5" system: More than 8 hours; Verso 8" system: More than 8 hours.


5. Cleaning

It is important to follow the cleaning procedures before using the equipment on human bodies; otherwise, patients and/or doctors may have the possibility of getting infection.

Wipe the pump unit with a damp cloth and a mild detergent, and keep it away from dust. If other detergent is used, choose one that will have no chemical effects on the surface of the plastic case of the pump unit.

 **Caution:** Do not immerse or soak pump unit.

Wipe down the mattress with warm water containing a mild detergent. The cover may also be cleaned by using sodium hypochlorite diluted in water. All parts should be air dried thoroughly before use.


 **CAUTION:** Do not use phenolic based product for cleaning.

 **CAUTION:** Keep mattresses out of direct exposure to sunlight, and make sure to dry the mattress after cleaning.

Washing Information Of Top Cover

Cover Material: Nylon/PU	    
Cover Material: Stretch	    

 **CAUTION:** Make sure to execute the cleaning procedure every time before a new patient lying on the air mattress.

 **Note:** More professional and detailed cleaning information can be found in the Cleaning & Disinfection Procedures published by APEX Medical. Please contact authorized distributor or qualified service center for information.

6. Storage

1. To store the mattress, lay the mattress out flat.
2. Roll from the head end towards the foot end with CPR valve open.


NOTE: Do not fold, crease or stack the mattresses, and keep mattresses out of direct exposure to sunlight.

7. Maintenance

7.1 General


1. Check main power cord and plug for abrasions or excessive wear.
2. Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stored together correctly.
3. Check airflow from the quick connector. The airflow should alternate between each connector every half-cycle time if it's in the alternating mode.
4. Check the air hoses if there is any kinks or breaks. For replacement, please contact our local service center (dealer).

7.2 Air Filter Replacement

1. Open the air filter plate located at the back of pump.
2. Suggest to check and replace air filter regularly if environment is dirty.

7.3 Fuse Replacement

1. Disconnect the plug from mains power when a fuse is to be replaced.
2. Remove the cover of the fuse holder by means of a screwdriver. (If necessary, un-clip the hanger from the pump unit)
3. Insert a new fuse of the correct rating in, and replace the cover of the fuse holder back. The fuse should be T1AL/250V type and VDE approved.

8. Expected Service Life

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

9. Trouble Shooting

If your questions can't be answered with above information, please contact our local service center directly. They might require a technician to rectify the problem.

Problem	Solution
Power is not ON	<ul style="list-style-type: none"> • Check if the plug is connected to mains. • Check if there is any blown fuse.
Low pressure indicator is ON	<ul style="list-style-type: none"> • Check if the connection between quick connector to pump unit is tightly secured. • Check if all tubing connections along the mattress are secured. • Check if the CPR valve is sealed. • Check if there is any air leakage on the air cells.
Patient is bottoming out	<ul style="list-style-type: none"> • Pressure setting might be inadequate for the patient, adjust pressure to higher level and wait for a few minutes.
Mattress form is loose	<ul style="list-style-type: none"> • Check if all the snap buttons or straps of mattress are all securely fastened. • Check if the mattress is fixed to the bed frame by straps.
No air produced from some air outlets of the air tube connector	<ul style="list-style-type: none"> • This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time.

10. Technical Specifications:

Item		Specification
Power Supply (Note: See rating label on the product)		AC230V 50Hz,0.05A (for 230V system)
Fuse Rating		T1A1, 250V
Cycle time		8 min for 60 Hz / 9.6 min for 50 Hz
Dimension (L x W x H)		25 x 13.5 x 9.5 (cm) / 9.8" x 5.3" x 3.7"
Weight		1.7 Kg
Environment	Temperature	Operation:10° C to 40° C (50° F to 104° F) Storage: -15° C to 50° C (5° F to 122° F) Shipping: -15° C to 70° C (5° F to 158° F)
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping:10 % to 90% non-condensing
	Atmospheric Pressure	700 hPa to 1060 hPa
Classification		Class II, Type BF, IP21(for 230V system) Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)
Mattress		Specification
Model		5" Overlay 8" Replacement
Dimension (L x W x H)		200 x 85/90 x 13 cm 200 x 85/90 x 20 cm 78.7" x 33.4"/35.4 x 5" 78.7" x 33.4"/35.4 x 8"
Weight		Width: 85cm, 4Kg Width: 85cm, 6Kg Width: 90cm, 4.5 Kg Width: 90cm, 7Kg
Max. Support Weight		160 Kg 180 Kg
Pressure Range		30-60 mmHg

NOTE:

1. Consult the distributor or EU representative for further technical documents.
2. These specifications are also applicable for other regions operating with the same power supply.
3. The earthing terminal is served as functional earth, not protective earth.(for UL only)
4. Mattress dimensions and weight is measured without foam base.
5. The manufacturer reserves the right to modify the specifications without prior notice.

Appendix A: EMC Information

Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
 Warning: <ol style="list-style-type: none"> 1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. 2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. 		

Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8kV contact ±15kV air		±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/output line		±2kV for power supply line ±1kV for input/output line	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		230V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the commended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF EM Fields IEC61000-4-3	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz 385-6000 MHz,	10 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz 385-6000 MHz,	10V/m	Recommended separation distance $d = \sqrt{P}$ 150kHz to 80MHz $d = 0.6\sqrt{P}$ 80MHz to 800MHz $d = 1.2\sqrt{P}$ 800 MHz to 2.7GHz

	<p>9-28V/m, 80% AM(1kHz) pulse mode and other modulation</p>	<p>9-28V/m, 80% AM(1kHz) pulse mode and other modulation</p>	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency ranged.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1: U_T is the a.c. mains voltage prior to the application of the test level</p> <p>NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people</p>			
<p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as repositioning or relocating the device.</p> <p>b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2.7 GHz $d = 1.2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2021/All rights reserved
Apex Medical (Kunshan) Corp.