

Upholstery Care

- **ALWAYS** ensure that you follow the cleaning instructions for **YOUR PARTICULAR UPHOLSTERY**. Your supplier will provide you with cleaning instructions. Please note that different fabrics will have different characteristics and cleaning instructions. Ask your supplier if you are in any doubt about the cleaning instructions. Repose cannot be held responsible for any damage caused by cleaning.
- Small dust and dirt particles can increase the wear and tear of the fabric and dull the surface appearance. Therefore, we recommend that you either brush your furniture using a soft brush or lightly vacuum. **NEVER** use abrasive tools or cleaning fluids.
- Your fabric should not be exposed to direct sunlight for long periods of time as this will cause your fabric to fade. If your furniture is positioned within prolonged light, ensure your curtains and blinds are drawn for protection.
- With certain fabrics such as Draylons and Chenilles pile flattening may occur, which is not a problem, but an inherent characteristic. This will appear as different shades and may give the impression of the fabric fading or creasing, which is caused by the pile of the fabric being moved in different directions. If this does occur, we recommend gently brushing in the opposite direction or vacuuming and leaving it to settle for a while. In most cases this will bring the pile back to its original position. However, you may experience some marking afterwards, which is unavoidable due to the nature of the fabric.
- Ensure your upholstery is protected from sharp objects or buckles as these may cause the upholstery to snag or become damaged. If the fabric snags do not pull the thread, carefully trim the thread with scissors.
- Avoid allowing pets onto your furniture as their claws may cause damage to the upholstery. All animals have a degree of oil in their coats which may transfer to your fabric if they go onto the chair. This will cause your fabric to prematurely wear.
- If you have fibre filled cushions regularly plump to avoid the fibres flattening.
- You may wish to protect certain areas of your chair by purchasing additional arm, head and seat covers from your supplier.

Warranty Upgrades

Certain chair usage conditions apply to these warranty upgrades. Check with your supplier to see whether you are eligible.

Extended Fabric Warranty / Accidental Damage Cover

You have 7 days after your delivery has taken place to decide whether you would like to purchase an additional warranty through your supplier. If purchased you will receive a small fabric treatment spray, Easy Care Cleaning Kit and a 5 year Fabric Protection and Accidental Damage Warranty.

Extended Structural Warranty

You have 7 days after your delivery has taken place to decide whether you would like to purchase an additional structural warranty. Ask your supplier about this 5 year Fabric Protection and Accidental Damage warranty which covers all internal and structural parts of your chair.

Warranty

Your Repose chair comes with a 2 year warranty on the mechanics, metal frame and electrical movement and handset. You also have a 1 year warranty on the wooden frame and upholstery.

Repose will repair or replace, at our discretion and to the original purchaser, any parts covered by this warranty and that on examination, are found to be defective in material and/or workmanship.

The warranty is for parts only and for transportation to and from the original delivery point.

Warranty Exclusions

- Fabric damage resulting from misuse, negligence, abuse, accident or incorrect cleaning.
- Unauthorised repairs, alterations and modifications including dismantling the chair (other than removing or fitting the back in line with these instructions).
- Damage to the chair from misuse, negligence, abuse or accident.
- Damage caused by common carrier abuse or mishandling and any special, indirect, incidental and/or consequential damages.

Warranty Claims

Please contact your supplier with your original purchase invoice and unique chair batch number and they will make the necessary arrangements with Repose.



Thank you for supporting British manufacturers
If you are in any doubt about any information or instructions
please contact your local supplier for further assistance.

Batch Number:

Recommended supplier:

Repose Furniture Ltd reserves the right to modify or change the specification without prior notification. For terms and conditions of sale and/or agreement please refer to your supplier.

www.reposefurniture.co.uk

Reposers

“Chairs designed with your personal comfort and wellbeing in mind”



Your Chair Getting Started

Your essential information on:

- Safety Instructions
- Chair Set Up
- Handset Instructions
- Upholstery Care
- Warranty and Contact Details



COMFORT WITHOUT COMPROMISE

Thank you for purchasing a Repose chair. Please take time to read this leaflet which contains essential instructions on your chair set up, handset usage and upholstery care as well as very important safety, warranty and contact details. However, if for any reason you are unsure of anything, please contact your supplier immediately for further information.

Important Safety Instructions

The following instructions are for your own safety, the safety of others and to ensure that the chair is used correctly and in line with the warranty conditions.

PLEASE ENSURE THAT:

- **YOU KEEP THESE INSTRUCTIONS SAFE** as you will need them for reference purposes.
- **YOU UNDERSTAND** the recommended user weight limit for your chair. It will be dependent on the Repose chair model that you have chosen. Please ensure that you understand these limits and **DO NOT** exceed them. This may damage the chair and **will invalidate your warranty.**
- **YOU CHECK** the chair and cables regularly to ensure that they are in a suitable working condition and that there is no damage. **DO NOT USE** the chair if it is not working properly or damaged.
- **YOU ALWAYS** use the base of the chair to move it around. **DO NOT** use the backrest to move the chair as this may cause damage to the framework and workings of the chair.
- **NO ONE** is sitting in the chair when it is being moved.
- **YOU ALWAYS** store your handset in either the pocket or on the loop provided. This will ensure that it will not be damaged by either falling on the floor or from anyone sitting on it.
- **NO** cables are trapped under the chair and that **YOU DO NOT** run over the cables with the castors or other equipment e.g. vacuum cleaner.
- Children **DO NOT** play with the handset and on, or with, the chair. To do so may cause personal injury.
- **NO ONE** stands or sits on the arms or leg rest as this may cause the chair to tip and/or damage the chair. To do so may cause personal injury.
- **NO ONE** stands on the seat as this may cause the chair to tip and/or damage the chair. To do so may cause personal injury.
- **YOU NEVER** exit the chair whilst the leg rest is elevated.
- **YOU ONLY USE** approved parts on your chair. These can be obtained through your supplier. **Failure to use approved parts will invalidate your warranty.**
- **YOU ALWAYS** get your supplier to carry out any repairs or fit new parts to the chair. **DO NOT** attempt to carry out your own repairs to the chair as this **will invalidate your warranty.**
- **YOU ALWAYS** disconnect your chair from the power source before checking any part of the chair.
- **YOU PROVIDE** your original purchase invoice and unique chair batch number when ordering parts or when making a claim under your warranty.

Delivering
Specialist Seating Solutions

Chair Location

- Position the chair in your chosen location, ensuring that the chair has enough space to recline fully and freely. There should be no objects behind the chair that will obstruct the chair movement and it should never touch the wall when it reclines. Such obstructions will cause damage to the chair and its mechanism and **will invalidate your warranty.**
- The leg rest should move freely with no objects to restrict this movement.
- Position the chair free from any trailing wires.
- Ensure that the chair is within easy reach of a power point to plug in the transformer (unless fitted with an Accupak) and the air pump.

Chair Set Up

- To position the back onto the base of the chair stand behind the chair and lay the back of the chair, face down and on top of the arms, ensuring that you align the brackets on the back rest to the upright supports on the base of the chair. Using slight pressure slot the back onto these supports. Depending on the mechanism you may find a small lever attached to both brackets. These should be pulled towards you to lock the back into position. Fit the back flap by using the Velcro attachments to complete the back set up.
- If you need to remove the back, release the locking bars (where applicable) and pull the back upwards and off.
- If you have an electric powered chair for safety reasons check that the extension cable, which is located at the rear of the chair, has not been damaged in transportation. On satisfaction connect the chair to your transformer.
- **ALWAYS** ensure that the transformer is placed in a flat position and is the correct way up.
- Once the transformer has been connected, plug the unit into the power source and switch the power on ready for use. **NB: The chair is safe when the power source is left on.**

Handset Instructions


Handsets come in varying designs depending on the type of mechanism inside the chair. The following instructions will aid easy use of your chair.

Single Motor (Dewert/Okini)

Button 1 – Back up, leg rest in, chair rises
Button 2 – Leg rest out, back down, chair lowers

Dual Motor (Dewert/Okini)

Button 1 – Back up
Button 2 – Chair rises, leg rest in
Button 3 – Back down
Button 4 – Chair lowers, leg rest out



Dual Motor (Okini)

Button 1 – Back down
Button 2 – Back up
Button 3 – Chair lowers, leg rest out
Button 4 – Chair rises, leg rest in
Button 5 – Chair lowers, leg rest out
Button 6 – Chair rises, leg rest in



Dual Motor (Stawert)

Button 1 – Chair lowers, leg rest out
Button 2 – Chair rises, leg rest in
Button 3 – Back down
Button 4 – Back up



Single Motor (Limoss)

Button 1 – Back up, leg rest in, chair rises
Button 2 – Leg rest out, back down, chair lowers



Dual Motor (Limoss)

Button 1 – Leg rest in, chair rises
Button 2 – Leg rest out, chair lowers
Button 3 – Back down
Button 4 – Back up




Single Motor (Okini)

Button 1 – Chair rises, leg rest in
Button 2 – Chair lowers, back down, leg rest out
Button 3 – Handset activated by pressing once. Press a second time the lights go off and handset locks. If you do not press a second time the handset will automatically go to lock down after a short period of time. This is for safety reasons.



Dual Motor (Okini)

Button 1 – Back up
Button 2 – Leg rest in, chair rises
Button 3 – Back down
Button 4 – Chair lowers, leg rest out
Button 5 – Handset activated by pressing once. Press a second time the lights go off and handset locks. If you do not press a second time the handset will automatically go to lock down after a short period of time. This is for safety reasons.



For heat or massage systems please see separate instruction leaflet.

Basic Technical Checks

ALWAYS disconnect your chair from the power source before checking any part of the chair. Then

- Check the wall socket is in working order by trying another appliance.
- Check the connection from the chair to the transformer is connected.
- Check to see if any of the cables are damaged or trapped. Damaged cables should be reported to your supplier.
- Check the electrical fuse and batteries (where applicable.) **DO NOT** open or attempt to open any electrical components other than to fit, as instructed, batteries or electrical fuses. **To do so may invalidate your warranty.**

www.reposefurniture.co.uk

COMFORT WITHOUT COMPROMISE