




**Rotunda Patient Transfer Aid**  
[www.enableaid.co.uk](http://www.enableaid.co.uk)

# How it works



The assistant wheels the Rotunda into position in front of the seated patient and ensures the platform turntable is **locked** and the wheels are facing the assistant.



The patient places their feet onto the turntable, they place their hands on the handle grip area and are assisted from a seated to a standing position.



The assistant unlocks and rotates the turntable smoothly left or right. There are two locking stations in either direction, once **locked** the patient is ready for transfer.



The assistant can now transfer the patient to a previously placed wheelchair, commode, chair or bed, the assistant then removes the Rotunda by tilting it back on to the wheels and moving it away.



**THE ROTUNDA IS A  
TRANSFER AID ONLY AND  
MUST NOT BE USED TO  
TRANSPORT PATIENTS**

# Accessories to Compliment your Rotunda

All accessories can be retro fitted and are suitable for all models.

## Knee Support

RAKS1



An easy clean, moulded Knee Support will provide extra comfort and aid weight bearing.

## Raising Blocks

76mm (3") RARB3 127mm (5") RARB5



Available in 76mm (3") or 127mm (5") heights, the Raising Blocks can be fitted together in any combination.

## Arm Supports

RAASH1



To aid weight bearing.

## Foot Guard

RAFG1



Foot Guard with instruction for use graphics.

Also available: Mat with Feet Placement Graphics RAMFP1. An easy clean mat with self adhesive backing.

## Cleaning/Maintenance Instructions

To clean, use a soft cloth and a solution of mild detergent and warm water. Do not use abrasive cleaners. The Rotunda needs no maintenance but should be inspected on a regular basis for any signs of wear and tear. The Rotunda is not designed for use in a wet room.

## One Year Warranty

Thank you for purchasing a Rotunda which has been designed manufactured and tested in the UK to the highest standards by Enable Aid Limited.

This Warranty covers defects in materials or workmanship for a period of one year from the date of original purchase.

During the Warranty period Enable Aid Limited will replace or repair, at no cost, any part of the Rotunda which proves to be defective due to manufacture defects or poor workmanship. If Enable Aid Limited is unable to provide a replacement part or undertake a suitable repair or if such repair is impracticable or cannot be made in a timely fashion Enable Aid Limited reserves the right to replace the whole unit.

This Warranty does not cover; Damage not resulting from defects in materials or workmanship, wear and tear, improper usage, misuse or abuse, negligence, accident or alteration.

This does not affect your statutory rights.

To make a claim under this warranty, in the first instance, please contact:


### **Enable Aid Limited**

Hawthorn Works  
33 Hollands Road  
Haverhill, Suffolk CB9 8PU

**Tel: 01440 702811** or

**Email: [sales@enableaid.co.uk](mailto:sales@enableaid.co.uk)**

Stating the nature of the problem and the date of purchase.



**We hope you enjoy your Rotunda and have many years trouble free usage**