

Tell us what you think

Your feedback is very important to us, so please take a moment and comment on our product and service.



Fast delivery on all UK orders

We deliver our products quickly and efficiently through courier or our specialist delivery team.



Head Office:

Copenhagen Road, Sutton Fields Industrial Estate, Kingston Upon Hull, East Yorkshire, HU7 0XQ.

E - sales@orchid-medicare.co.uk

- T 01482 820333
- F 01482 824617



Related Accessories

The tilt -in space shower chair has an extensive range of optional accessories designed to assist disabled people and carers in their everyday needs.

| Commode Facility | Detachable Side Panels | Head Rest | Mesh Covered Neck Rolls |
|---|------------------------|---------------------|----------------------------|
| Mesh Covered Foam Side Support Rolls 4 x12' (Pair) | Protective Mesh Cover | Adjustable footrest | Overlay Panel for Aperture |
| Chest Straps | Full Body Harness | Ankle Straps | Knee Straps |
| | | | |





ORCHID ΤΙΔΤ-ΙΝ S P A C E **SHOWER CHAIR** www.orchid-medicare.co.uk

Designed to assist MADE IN THE UK







The Orchid tilt in space shower chair has been designed to provide comfort and safety together with easy access when washing and for personal care.

The shower chair has optional extra such as fully adjustable and removable footrest, thoracic supports and hip supports. The shower chair is available in both adult and paediatric sizes and the option of either manual and electric tilting actions both of which provide a tilting angle between 75° and 35°.

Both of the chairs are available with or without an aperture. The frame is made of stainless steel so will not rust and there is the added protection of powder coating with BacAway an anti-microbial anti-fungal agent preventing infection.



OPTIONAL ACCESSORIES

■ Velcro Rolls / Support Rolls



Head Rest



TILT-IN SPACE SHOWER CHAIR **Designed to assist**

Detachable side panels provide additional security for the client. Easy to use tilt action Adjustable control with smooth and safety belt. reliable operation. Base transporter stainless steel frame with a powder coat finish. A quick release mesh Dual gas for ease of laundry or spring replacement. assembly. Heavy duty castor with easy to use Adjustable foot rest brake, made from for individual comfort. high grade materials.

See our wide range of accessories and replacement products





Easy control reclining action gives the patient maximum comfort



■ The top of the mesh can be released

rear to allow optimum access for

hair washing.

by undoing the velcro straps at the









TECHNICAL DATA

No specialist cleaning required due to our BacAway anti-bacterial powder coating, just use a proprietary cleaning solution.

Maximum Overall Length Reclined (with footrest) = Adult 1480mm / Child 1200mm.

Maximum Overall Length Upright (with footrest) = Adult 1010mm / Child 920mm.

Overall Width = 670mm.

Overall Height (in upright position) = Adult 1360mm / Child 1160mm.

Overall Height (reclined)


= Adult 910mm / Child 780mm.

Small Chair Weight = 22kg

Large Chair Weight = 23kg

Orchid Medicare - Tilt-in Space Shower Chair Standard Sizing Matrix Note Dimensions A-D are overall widths

| Dimension | Adult frame sizes | Child frame sizes |
|------------------------------|-------------------|-------------------|
| A - (Back Width) | 573mm / 221/2" | 500mm / 19¾" |
| B - (Hip) | 510mm / 20" | 455mm / 18" |
| C - (Knee) | 540mm / 21¼" | 495mm / 19½" |
| D - (Foot) | 405mm / 16" | 405mm / 16" |
| E - (Back Length) | 812mm / 32" | 612mm / 24" |
| F - (Thigh) | 405mm / 16" | 355mm / 14" |
| G - (Calf) | 335mm / 131⁄4" | 260mm / 10¼" |
| H - (Overall Width) | 670mm / 26½" | 670mm / 26½" |
| I - (Seat Ht) | 580mm / 22¾" | 580mm / 22¾" |
| J - (Overall Length Upright) | 1010mm / 39¾" | 920mm / 36¼" |
| K - (Overall Ht) | 1360mm / 53½" | 1160mm / 45¾" |
| Overall Length Reclined | 1480mm / 58¼" | 1200mm / 471/4" |



CONSTRUCTION

- Tubular stainless steel material ensuring robust contruction.
- Our Anti bacterial powder coating which does not degrade effectiveness over time.
- Castor tested to 100 hours salt spray test.
- Upholstery is hard wearing easy clean and meets current legislations regarding fire retardant upholstery.