













Slowing brake

The "Sparrow" is a lightweight, easy to use wheelchair which offers outstanding value for a durable product (Karma provide two years warranty on the frame for added peace of mind). The transit version weighs just 11.7kg (25lbs) complete and easily folds into a small package for transport/storage.

Features

- Lightweight from 11.7kg
- Maximum user weight 100kgs
- Folding frame/back
- Puncture proof tyres
- Slowing brakes

Standard specification

- Standard 17"seat width
- 8"solid castors
- Fixed 14" Mag Wheels/Quick release 24" rear wheels
- Fixed armrest
- Fixed footrest
- Black flame retardant upholstery





tions	ite Model	ems Width	F/R wheels		St.	<u> </u>			<u></u>	<u></u>		<u>b.</u>		Max. User Weight
Specificat	SM-150 Q24	17"	8"-24"	42	43	68	29	23	100	48.5	43.5	92	13,15	100kg
	SM-150 F14	17"	8"-14"	42	43	66	27	23	87	48.5	43.5	92	11.79	100kg
	KARMA reserves the right to modify information herein without further notice. •												<u>+</u> 1cn	n <u>+</u> 0.5kg

KARMA Mobility Limited

104-40537 Rev.02 OCT.2008